

Product Datasheet

Mouse Secreted Frizzled Related Protein 1 ELISA Kit (orb1212482)

Description	Mouse Secreted Frizzled Related Protein 1, SFRP1 ELISA Kit
Reactivity	Mouse
Range	10-2000ng/mL
Tested Applications	ELISA
Target	SFRP1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	35,413
Uniprot ID	Q8C4U3
Solubility (25°C)	2400ng/mL
Sensitivity	5.29ng/mL
Entrez	20377
Expiration Date	Please enquire.