

Product Datasheet

Mouse Secreted Frizzled Related Protein 2 ELISA Kit (orb1212483)

Description	Mouse Secreted Frizzled Related Protein 2, SFRP2 ELISA Kit
Reactivity	Mouse
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	SFRP2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	33,469
Uniprot ID	P97299
Solubility (25°C)	48ng/ml
Sensitivity	0.054ng/ml
Entrez	20319
Expiration Date	Please enquire.