

Product Datasheet

Mouse Vascular Endothelial Cell Growth Factor C ELISA Kit (orb1212514)

Description	Mouse Vascular Endothelial Cell Growth Factor C, VEGF-C ELISA Kit
Reactivity	Mouse
Range	0.5-200ng/L
Tested Applications	ELISA
Target	VEGFC
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	46,471
Uniprot ID	P97953
Solubility (25°C)	240ng/L
Sensitivity	0.31ng/L
Entrez	22341
Expiration Date	Please enquire.