

Product Datasheet

Mouse Insulin-like Growth Factor 2 ELISA Kit (orb1212523)

Description	Mouse Insulin-like Growth Factor 2, IGF-2 ELISA Kit
Reactivity	Mouse
Range	5-1500ng/ml
Tested Applications	ELISA
Target	IGF2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	20,030
Uniprot ID	P09535
Solubility (25°C)	1600ng/ml
Sensitivity	2.28ng/ml
Entrez	16002
Expiration Date	Please enquire.