

Product Datasheet

Mouse Complement Factor H ELISA Kit (orb1212534)

Description	Mouse Complement Factor H, CFH ELISA Kit
Reactivity	Mouse
Range	20-4500ng/mL
Tested Applications	ELISA
Target	CFH
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	139,138
Uniprot ID	P06909
Solubility (25°C)	4800ng/mL
Sensitivity	9.76ng/mL
Entrez	12628
Expiration Date	Please enquire.