

Product Datasheet

Mouse Hyaluronan Synthase 2 ELISA Kit (orb1212537)

Description	Mouse Hyaluronan Synthase 2, HAS2 ELISA Kit
Reactivity	Mouse
Range	15-3000ng/L
Tested Applications	ELISA
Target	HAS2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	63,510
Uniprot ID	P70312
Solubility (25°C)	3200ng/L
Sensitivity	6.94ng/L
Entrez	15117
Expiration Date	Please enquire.