

Product Datasheet

Mouse Toll-like Receptor 2 ELISA Kit (orb1212548)

Description	Mouse Toll-like Receptor 2, TLR2 ELISA Kit
Reactivity	Mouse
Range	0.05-15ng/ml
Tested Applications	ELISA
Target	TLR2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	89,449
Uniprot ID	Q9QUN7
Solubility (25°C)	16ng/ml
Sensitivity	0.01ng/ml
Entrez	24088
Expiration Date	Please enquire.