

Product Datasheet

Mouse Substance P Receptor ELISA Kit (orb1212552)

Description	Mouse Substance P Receptor, SP-R ELISA Kit
Reactivity	Mouse
Range	3-900ng/L
Tested Applications	ELISA
Target	TACR1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	46,322
Uniprot ID	P30548
Solubility (25°C)	960ng/L
Sensitivity	1.53ng/L
Entrez	21336
Expiration Date	Please enquire.