

Product Datasheet

Mouse Insulin-like Growth Factor Binding Protein 2 ELISA Kit (orb1212574)

Description	Mouse Insulin-like Growth Factor Binding Protein 2, IGFBP-2 ELISA Kit
Reactivity	Mouse
Range	0.2-60ng/mL
Tested Applications	ELISA
Target	IGFBP2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	32,847
Uniprot ID	P47877
Solubility (25°C)	64ng/mL
Sensitivity	0.131ng/mL
Entrez	16008
Expiration Date	Please enquire.