

Product Datasheet

Mouse Thrombopoietin ELISA Kit (orb1212587)

Description	Mouse Thrombopoietin, TPO ELISA Kit
Reactivity	Mouse
Range	3-900ng/L
Tested Applications	ELISA
Target	THPO
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	37,836
Uniprot ID	P40226
Solubility (25°C)	960ng/L
Sensitivity	1.54ng/L
Entrez	21832
Expiration Date	Please enquire.