

Product Datasheet

Mouse Interleukin 5 ELISA Kit (orb1212601)

Description	Mouse Interleukin 5, IL-5 ELISA Kit
Reactivity	Mouse
Range	1-400ng/L
Tested Applications	ELISA
Target	IL5
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	15,410
Uniprot ID	P04401
Solubility (25°C)	480ng/L
Sensitivity	0.51ng/L
Entrez	16191
Expiration Date	Please enquire.