

## Product Datasheet

# Mouse Insulin-like Growth Factor 1 ELISA Kit (orb1212607)

<b>Description</b>	Mouse Insulin-like Growth Factor 1, IGF-1 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	0.05-10ng/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	IGF1
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	17,093
<b>Uniprot ID</b>	<b>P05017</b>
<b>Solubility (25°C)</b>	12.8ng/ml
<b>Sensitivity</b>	0.01ng/ml
<b>Entrez</b>	<b>16000</b>
<b>Expiration Date</b>	Please enquire.