

Product Datasheet

Mouse Interleukin 13 ELISA Kit (orb1212614)

Description	Mouse Interleukin 13, IL-13 ELISA Kit
Reactivity	Mouse
Range	1-400ng/L
Tested Applications	ELISA
Target	IL13
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	14,108
Uniprot ID	P20109
Solubility (25°C)	400ng/L
Sensitivity	0.52ng/L
Entrez	16163
Expiration Date	Please enquire.