

Product Datasheet

Human 40S Ribosomal Protein S23 ELISA Kit (orb1212637)

Description	Human 40S Ribosomal Protein S23, RPS23 ELISA Kit
Reactivity	Human
Range	7-1500ng/L
Tested Applications	ELISA
Target	RPS23
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	15,808
Uniprot ID	P62266
Solubility (25°C)	1600ng/L
Sensitivity	3.55ng/L
Entrez	6228
Expiration Date	Please enquire.