

## Product Datasheet

# Human Eukaryotic Translation Initiation Factor 2-alpha Kinase 3 ELISA Kit (orb1212653)

<b>Description</b>	Human Eukaryotic Translation Initiation Factor 2-alpha Kinase 3, EIF2AK3 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.05-15ng/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	EIF2AK3
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	125,216
<b>Uniprot ID</b>	<a href="#">Q9NZJ5</a>
<b>Solubility (25°C)</b>	16ng/ml
<b>Sensitivity</b>	0.019ng/ml
<b>Entrez</b>	<a href="#">9451</a>
<b>Expiration Date</b>	Please enquire.