

Product Datasheet

Human Orexin ELISA Kit (orb1212655)

Description	Human Orexin, HCRT ELISA Kit
Reactivity	Human
Range	37.5-1200pg/ml
Tested Applications	ELISA
Target	HCRT
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	13,363
Uniprot ID	O43612
Solubility (25°C)	2400pg/ml
Sensitivity	15.42pg/ml
Entrez	3060
Expiration Date	Please enquire.