

## Product Datasheet

# Human Microtubule-associated Proteins 1A ELISA Kit (orb1212684)

<b>Description</b>	Human Microtubule-associated Proteins 1A, 1B Light Chain 3A, MAP1LC3A ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	20-6000ng/L
<b>Tested Applications</b>	ELISA
<b>Target</b>	MAP1LC3A
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	14,272
<b>Uniprot ID</b>	<a href="#">Q9H492</a>
<b>Solubility (25°C)</b>	6400ng/L
<b>Sensitivity</b>	11.63ng/L
<b>Entrez</b>	<a href="#">84557</a>
<b>Expiration Date</b>	Please enquire.