

Product Datasheet

Human Tumor Necrosis Factor Receptor Superfamily Member 1B ELISA Kit (orb1212686)

Description	Human Tumor Necrosis Factor Receptor Superfamily Member 1B, TNFRSF1B ELISA Kit
Reactivity	Human
Range	20-3800ng/L
Tested Applications	ELISA
Target	TNFRSF1B
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	48,291
Uniprot ID	P20333
Solubility (25°C)	4000ng/L
Sensitivity	10.98ng/L
Entrez	7133
Expiration Date	Please enquire.