

Product Datasheet

Human 5-Hydroxytryptamine Receptor 2A ELISA Kit (orb1212712)

Description	Human 5-Hydroxytryptamine Receptor 2A, HTR2A ELISA Kit
Reactivity	Human
Range	0.05-10ng/ml
Tested Applications	ELISA
Target	HTR2A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	52,603
Uniprot ID	P28223
Solubility (25°C)	12.8ng/ml
Sensitivity	0.019ng/ml
Entrez	3356
Expiration Date	Please enquire.