

Product Datasheet

Human 5-Hydroxytryptamine Receptor 3A ELISA Kit (orb1212714)

Description	Human 5-Hydroxytryptamine Receptor 3A, HTR3A ELISA Kit
Reactivity	Human
Range	20-3800ng/L
Tested Applications	ELISA
Target	HTR3A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	55,280
Uniprot ID	P46098
Solubility (25°C)	4000ng/L
Sensitivity	11.55ng/L
Entrez	3359
Expiration Date	Please enquire.