

Product Datasheet

Human Insulin Like Growth Factor Binding Protein 4 ELISA Kit (orb1212726)

Description	Human Insulin Like Growth Factor Binding Protein 4, IGFBP4 ELISA Kit
Reactivity	Human
Range	1-400ng/ml
Tested Applications	ELISA
Target	IGFBP4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	27,934
Uniprot ID	P22692
Solubility (25°C)	480ng/ml
Sensitivity	0.53ng/ml
Entrez	3487
Expiration Date	Please enquire.