

Product Datasheet

Human Tumor Necrosis Factor Receptor Superfamily Member 17 ELISA Kit (orb1212742)

Description	Human Tumor Necrosis Factor Receptor Superfamily Member 17, TNFRSF17 ELISA Kit
Reactivity	Human
Range	3-900ng/L
Tested Applications	ELISA
Target	TNFRSF17
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	20,165
Uniprot ID	Q02223
Solubility (25°C)	960ng/L
Sensitivity	1.62ng/L
Entrez	608
Expiration Date	Please enquire.