

Product Datasheet

Human Cytochrome P450 2C19 ELISA Kit (orb1212748)

Description	Human Cytochrome P450 2C19, CYP2C19 ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	CYP2C19
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	55,931
Uniprot ID	P33261
Solubility (25°C)	48ng/ml
Sensitivity	0.053ng/ml
Entrez	1557
Expiration Date	Please enquire.