

Product Datasheet

Human Calcyclin-binding Protein ELISA Kit (orb1212754)

Description	Human Calcyclin-binding Protein, CACYBP ELISA Kit
Reactivity	Human
Range	20-3800ng/L
Tested Applications	ELISA
Target	CACYBP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	26,210
Uniprot ID	Q9HB71
Solubility (25°C)	4000ng/L
Sensitivity	9.17ng/L
Entrez	27101
Expiration Date	Please enquire.