

Product Datasheet

Human V-Set and Immunoglobulin Domain-containing Protein 4 ELISA Kit (orb1212803)

Description	Human V-Set and Immunoglobulin Domain-containing Protein 4, VSIG4 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	VSIG4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	43,987
Uniprot ID	Q9Y279
Solubility (25°C)	3200ng/L
Sensitivity	7.69ng/L
Entrez	11326
Expiration Date	Please enquire.