

Product Datasheet

Human Tubulin Beta-3 Chain ELISA Kit (orb1212815)

Description	Human Tubulin Beta-3 Chain, TUBB3 ELISA Kit
Reactivity	Human
Range	0.3-90ng/ml
Tested Applications	ELISA
Target	TUBB3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	50,433
Uniprot ID	Q13509
Solubility (25°C)	96ng/ml
Sensitivity	0.146ng/ml
Entrez	10381
Expiration Date	Please enquire.