

## Product Datasheet

# Human Thioredoxin-dependent Peroxide Reductase ELISA Kit (orb1212842)

<b>Description</b>	Human Thioredoxin-dependent Peroxide Reductase, Mitochondrial, PRDX3 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	2-700ng/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	PRDX3
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	27,693
<b>Uniprot ID</b>	<b>P30048</b>
<b>Solubility (25°C)</b>	800ng/mL
<b>Sensitivity</b>	0.91ng/mL
<b>Entrez</b>	<b>10935</b>
<b>Expiration Date</b>	Please enquire.