

Product Datasheet

Human Sortilin-related Receptor ELISA Kit (orb1212861)

Description	Human Sortilin-related Receptor, SORL1 ELISA Kit
Reactivity	Human
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	SORL1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	248,426
Uniprot ID	Q92673
Solubility (25°C)	24ng/mL
Sensitivity	0.029ng/mL
Entrez	6653
Expiration Date	Please enquire.