

Product Datasheet

Human Prokineticin-2 ELISA Kit (orb1212944)

Description	Human Prokineticin-2, PROK2 ELISA Kit
Reactivity	Human
Range	3-900ng/L
Tested Applications	ELISA
Target	PROK2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	14,314
Uniprot ID	Q9HC23
Solubility (25°C)	960ng/L
Sensitivity	1.53ng/L
Entrez	60675
Expiration Date	Please enquire.