

Product Datasheet

Human Profilin 1 ELISA Kit (orb1212947)

Description Human Profilin 1 ELISA Kit

Reactivity Human

Range 20-4500ng/L

Tested Applications ELISA

Target PFN1

Storage 2-8°C

Note For research use only

Application notes Sandwich

Sample Types Serum, plasma, cell culture supernates

MW 15,054

Uniprot ID [P07737](#)

Solubility (25°C) 4800ng/L

Sensitivity 9.33ng/L

Entrez [5216](#)

Expiration Date Please enquire.