

## Product Datasheet

# Human Peptidoglycan Recognition Protein 1 ELISA Kit (orb1212971)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Human Peptidoglycan Recognition Protein 1, PGLYRP1 ELISA Kit |
| <b>Reactivity</b>          | Human  |
| <b>Range</b>               | 10-2000ng/L  |
| <b>Tested Applications</b> | ELISA  |
| <b>Target</b>              | PGLYRP1  |
| <b>Storage</b>             | 2-8°C  |
| <b>Note</b>                | For research use only  |
| <b>Application notes</b>   | Sandwich   |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates                       |
| <b>MW</b>                  | 21,731   |
| <b>Uniprot ID</b>          | <a href="#">O75594</a>                                       |
| <b>Solubility (25°C)</b>   | 2400ng/L   |
| <b>Sensitivity</b>         | 4.28ng/L   |
| <b>Entrez</b>              | <a href="#">8993</a>   |
| <b>Expiration Date</b>     | Please enquire.  |