

Product Datasheet

Human Pappalysin-2 ELISA Kit (orb1212974)

Description	Human Pappalysin-2, PAPPA2 ELISA Kit
Reactivity	Human
Range	0.5-100ng/mL
Tested Applications	ELISA
Target	PAPPA2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	198,539
Uniprot ID	Q9BXP8
Solubility (25°C)	128ng/mL
Sensitivity	0.215ng/mL
Entrez	60676
Expiration Date	Please enquire.