

Product Datasheet

Human Myelin-associated Glycoprotein ELISA Kit (orb1213016)

Description	Human Myelin-associated Glycoprotein, MAG ELISA Kit
Reactivity	Human
Range	0.3-90ng/mL
Tested Applications	ELISA
Target	MAG
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	69,069
Uniprot ID	P20916
Solubility (25°C)	96ng/mL
Sensitivity	0.146ng/mL
Entrez	4099
Expiration Date	Please enquire.