

Product Datasheet

Human Multidrug Resistance-associated Protein 1 ELISA Kit (orb1213019)

| | |
|----------------------------|------------------------------------------------------------------|
| Description | Human Multidrug Resistance-associated Protein 1, ABCC1 ELISA Kit |
| Reactivity | Human |
| Range | 10-2000ng/mL |
| Tested Applications | ELISA |
| Target | ABCC1 |
| Storage | 2-8°C |
| Note | For research use only |
| Application notes | Sandwich |
| Sample Types | Serum, plasma, cell culture supernates |
| MW | 171,591 |
| Uniprot ID | P33527 |
| Solubility (25°C) | 2400ng/mL |
| Sensitivity | 4.34ng/mL |
| Entrez | 4363 |
| Expiration Date | Please enquire. |