

Product Datasheet

Human Microfibrillar-associated Protein 5 ELISA Kit (orb1213029)

Description	Human Microfibrillar-associated Protein 5, MFAP5 ELISA Kit
Reactivity	Human
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	MFAP5
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	19,612
Uniprot ID	Q13361
Solubility (25°C)	24ng/mL
Sensitivity	0.024ng/mL
Entrez	8076
Expiration Date	Please enquire.