

Product Datasheet

Human Interleukin-18 Receptor Accessory Protein ELISA Kit (orb1213087)

Description	Human Interleukin-18 Receptor Accessory Protein, IL18RAP ELISA Kit
Reactivity	Human
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	IL18RAP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	68,310
Uniprot ID	095256
Solubility (25°C)	25.6ng/mL
Sensitivity	0.028ng/mL
Entrez	8807
Expiration Date	Please enquire.