

Product Datasheet

Human Interleukin-18 Receptor 1 ELISA Kit (orb1213088)

Description	Human Interleukin-18 Receptor 1, IL18R1 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	IL18R1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	62,304
Uniprot ID	Q13478
Solubility (25°C)	25.6ng/ml
Sensitivity	0.026ng/ml
Entrez	8809
Expiration Date	Please enquire.