

Product Datasheet

Human Interleukin-17 Receptor C ELISA Kit (orb1213090)

Description	Human Interleukin-17 Receptor C, IL17RC ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	IL17RC
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	86,240
Uniprot ID	Q8NAC3
Solubility (25°C)	24ng/ml
Sensitivity	0.027ng/ml
Entrez	84818
Expiration Date	Please enquire.