

## Product Datasheet

# Human Interleukin-17 Receptor B ELISA Kit (orb1213091)

<b>Description</b>	Human Interleukin-17 Receptor B, IL17RB ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.3-90ng/L
<b>Tested Applications</b>	ELISA
<b>Target</b>	IL17RB
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	55,885
<b>Uniprot ID</b>	<a href="#">Q9NRM6</a>
<b>Solubility (25°C)</b>	96ng/L
<b>Sensitivity</b>	0.17ng/L
<b>Entrez</b>	<a href="#">55540</a>
<b>Expiration Date</b>	Please enquire.