

## Product Datasheet

# Human Insulin-like Growth Factor 2 mRNA-binding Protein 3 ELISA Kit (orb1213119)

<b>Description</b>	Human Insulin-like Growth Factor 2 mRNA-binding Protein 3, IGF2BP3 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	20-4500ng/L
<b>Tested Applications</b>	ELISA
<b>Target</b>	IGF2BP3
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	63,705
<b>Uniprot ID</b>	<a href="#">O00425</a>
<b>Solubility (25°C)</b>	4800ng/L
<b>Sensitivity</b>	9.87ng/L
<b>Entrez</b>	<a href="#">10643</a>
<b>Expiration Date</b>	Please enquire.