

Product Datasheet

Human Homeobox Protein Nanog ELISA Kit (orb1213136)

Description	Human Homeobox Protein Nanog, NANOG ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	NANOG
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	34,620
Uniprot ID	Q9H9S0
Solubility (25°C)	3200ng/L
Sensitivity	7.44ng/L
Entrez	79923
Expiration Date	Please enquire.