

Product Datasheet

Human Hepatitis A Virus Cellular Receptor 2 ELISA Kit (orb1213155)

Description	Human Hepatitis A Virus Cellular Receptor 2, HAVCR2 ELISA Kit
Reactivity	Human
Range	20-6000ng/L
Tested Applications	ELISA
Target	HAVCR2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	33,394
Uniprot ID	Q8TDQ0
Solubility (25°C)	6400ng/L
Sensitivity	9.95ng/L
Entrez	84868
Expiration Date	Please enquire.