

Product Datasheet

Human Heat-Stable Enterotoxin Receptor ELISA Kit (orb1213159)

Description	Human Heat-Stable Enterotoxin Receptor, GUCY2C ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	GUCY2C
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	123,403
Uniprot ID	P25092
Solubility (25°C)	3200ng/L
Sensitivity	7.97ng/L
Entrez	2984
Expiration Date	Please enquire.