

Product Datasheet

Human Glutamate Carboxypeptidase 2 ELISA Kit (orb1213179)

Description	Human Glutamate Carboxypeptidase 2, FOLH1 ELISA Kit
Reactivity	Human
Range	2-600ng/mL
Tested Applications	ELISA
Target	FOLH1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	84,331
Uniprot ID	Q04609
Solubility (25°C)	640ng/mL
Sensitivity	1.22ng/mL
Entrez	2346
Expiration Date	Please enquire.