

## Product Datasheet

### Human Frataxin ELISA Kit (orb1213195)

<b>Description</b>	Human Frataxin, Mitochondrial, FXN ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	15-3000ng/L
<b>Tested Applications</b>	ELISA
<b>Target</b>	FXN
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	23,135
<b>Uniprot ID</b>	<a href="#">Q16595</a>
<b>Solubility (25°C)</b>	3200ng/L
<b>Sensitivity</b>	7.81ng/L
<b>Entrez</b>	<a href="#">2395</a>
<b>Expiration Date</b>	Please enquire.