

Product Datasheet

Human Ets Domain-containing Protein Elk-3 ELISA Kit (orb1213218)

Description	Human Ets Domain-containing Protein Elk-3, ELK3 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	ELK3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	44,240
Uniprot ID	P41970
Solubility (25°C)	3200ng/L
Sensitivity	7.07ng/L
Entrez	2004
Expiration Date	Please enquire.