

Product Datasheet

Human Estrogen-related Receptor Gamma ELISA Kit (orb1213220)

Description	Human Estrogen-related Receptor Gamma, ESRRG ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	ESRRG
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	51,306
Uniprot ID	P62508
Solubility (25°C)	4800ng/L
Sensitivity	9.96ng/L
Entrez	2104
Expiration Date	Please enquire.