

Product Datasheet

Human Endothelin-converting Enzyme 2 ELISA Kit (orb1213230)

Description	Human Endothelin-converting Enzyme 2, ECE2 ELISA Kit
Reactivity	Human
Range	20-6000ng/L
Tested Applications	ELISA
Target	ECE2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	91,211
Uniprot ID	P0DPP6
Solubility (25°C)	6400ng/L
Sensitivity	9.74ng/L
Entrez	1105995839718
Expiration Date	Please enquire.