

## Product Datasheet

# Human DNA-binding Protein Inhibitor Id-2 ELISA Kit (orb1213242)

<b>Description</b>	Human DNA-binding Protein Inhibitor Id-2, ID2 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	15-3000ng/L
<b>Tested Applications</b>	ELISA
<b>Target</b>	ID2
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	14,917
<b>Uniprot ID</b>	<a href="#">Q02363</a>
<b>Solubility (25°C)</b>	3200ng/L
<b>Sensitivity</b>	9.48ng/L
<b>Entrez</b>	<a href="#">3398</a>
<b>Expiration Date</b>	Please enquire.