

Product Datasheet

Human C-X-C Motif Chemokine 3 ELISA Kit (orb1213275)

Description	Human C-X-C Motif Chemokine 3, CXCL3 ELISA Kit
Reactivity	Human
Range	5-1000ng/L
Tested Applications	ELISA
Target	CXCL3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	11,342
Uniprot ID	P19876
Solubility (25°C)	1280ng/L
Sensitivity	2.29ng/L
Entrez	2921
Expiration Date	Please enquire.